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REF 054 800

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D05480012022-041005 ML I Rev.01

DE BORT TalofX

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EN BORT TalofX

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is an adjustable ankle brace for the stabilisation of the upper and lower ankle joint on one plane. The brace consists of elastic material and stabilisation elements made of rigid material.

Indications

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory irritations, rheumatism, arthritis, tendinitis, ligament injuries, supination/pronation trauma.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy

Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material

Fitting instructions

1 Release the 3 fastening straps and loosen them. Slip the brace over the toes and put it on like a sock.

2 Then close the circular strap on the lower leg. Ensure that the brace doesn't twist on the foot and lower leg.

3 Then, close the stabilisation belt with the yellow marking. For this, guide it diagonally over the metatarsal externally through the free flap on the inner ankle towards the rear to the externally attached guide loop. Close the belt loosely there.

4 Then, close the stabilisation belt with the orange marking. For this, guide it over the metatarsal internally to the guide loop and close it there loosely. Then, tighten the harness closed in Step 3, then the belt strap last closed.

5 It is quite possible that the belt straps again need to be tightened after a few minutes have lapsed.

Removal

To remove the BORT TalofX, open the belts in the reverse order.

Tip: Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Information: The plastic stabilisation elements positioned at the side of the brace are not thermoplastically remouldable/mouldable.

Material composition

Polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer
 Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 04.2022

Medical device | Single patient – multiple use

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